

Superlite Electric Foldable wheelchair

Product Code: 9346376063944

Features:

- Carbon Fibre Frame
- Ultra Lightweight: 17KG (excluding battery) 20kg with battery
- Portable: folding size 86x41x71cm
- Easy to transport in a car boot
- 2 Brushless 180W Motors
- Battery 12Ah lithium
- > 25 km driving range on single charge
- Intelligent electromagnetic brake system
- Easy to fold and unfold





> 353710

Length (cm)	103	Backrest Height (cm)	14	Front Wheel Diameter (cm)	17.5cm
Armrest Width (cm)	45	Seat Depth (cm)	43	Rear Wheel Diameter (cm)	19.5cm
Height (cm)	94	Seat Width (cm)	42	Max Load (kg)	120kg
Folded size (cm)	41*86 *71	Seat Height from the Ground(cm)	50	Net Weight with- out battery (kg) Battery (kg)	17kg 3kg