
Healthcare furniture, enhancing lives.



EN7 Series Bed







The EN7 series bed is a variant of the EN9, manufactured in Australia from imported components. This bed is versatile, as is suitable for low, high, and bariatric care; featuring a low height of 100mm and a high height of 725mm.

The EN7 is available in single and king single sizes and is standard in an APO Grey frame and Natural Oak Laminate Head and Footboards.

The EN7 enables care for frail, elderly or dementia patients, through acute and bariatric settings, this bed has the needs of the residents and the carers covered. One bed, endless possibilities!

FEATURES

- o Extremely Low Minimum Position
- o Auto-Regression Backrest
- o 75mm Castors
- APO Grey Powdercoated Frame



Chair Position at the Push of a Button



Bluetooth Connectivity ***handset must be ordered separately***



Slender Arch Head & Footboard in Natural Oak



Vertical Lift

SPECIFICATIONS

EN7	Single	King Single
Overall Dimensions (mm)	2115 x 920	2115 x 1090
Sleeping Surface (mm)	2000 x 900	2000 x 1070
SWL (kg)	200	200
High/low Position (mm)	725/100	725/100
Back Angle (degrees)	62	62
Knee Angle (degrees)	35	35
Reverse Angle (degrees)	15	15
Trendelenburg Angle (degrees)	15	15
Product Codes	EN7S-00000	EN7KS-00000

WARRANTY: Frame - 10 years, Electrics - 5 years & Castors 1 year.



T: 1800 282 001

E: sales@caremed.healthcare W: www.caremed.healthcare



OPTIONAL FEATURES

Added at time of manufacture

- 300kg Handset 0
- Battery Back Up 0
- 0 Dual Breaking
- Head & Footboard Options



Light





Attendant Control



CPR Backrest



OPTIONAL HEAD & FOOTBOARDS

Handset

Custom options available

Laminate Slender Arch Style



Natural Oak

Upholstered Jersey Headboard & Plain Footboard



Ash Ocean

Timber Eltham Style



Burnt Walnut Beech Saddle

ACCESSORIES

- Mattresses
- Reading Light
- **Emergency Power Supply**





400HR-EN7



CGH Combro Grab Handle



ASHP Self Help Pole & Triangle



AIVP I.V Pole



EN9-TBR T-Bar Wall Protection



EN-SE-HH StandEzy Bed Pole



EN9-CLPDB Bedstick



C.EXTEN9000 Clip on Extension

T: 1800 282 001

E: sales@caremed.healthcare W: www.caremed.healthcare





Can't find what your looking for? Contact us to disucss your custom options.



Pick your colour



Choose your width



Measure required length



Create your style



We build for legends



with us

Distributed by:

