

Peak B-941 Electric Foldable Bed - Single

Product Code: 9346376013734

Features:

- ▶ Electric Folding, one person set up
- ▶ Compatible with BH989 Clamp on Bed Rail and self help pole
- ▶ 4 Locking Castor Wheels
- ▶ DC TiMotion Motor System
- ▶ Side rail and battery backup as option
- ▶ Steel Frame with Special Anti-Bacterial Coating



Product Specifications

Max Bed Height	70cm
Min Bed Height	30cm
Length	231cm
Width	100cm
Trendelenburg & Rev	-10° - +10°
Load Capacity	220kg
Product Weight	92kg

Product Specifications

Surface Length	198cm
Surface Width	86cm
Backrest Adjustment	0 - 73°
Knee- Break Adjustment	0 - 35°

ARTG Listing No.

- ▶ 196516

Peak B-941K Electric Foldable Bed - King Single bed

Product Code: 9346376014366

Features:

- Electric Folding, one person set up
- Compatible with BH989 Clamp on Bed Rail and self help pole
- 4 Locking Castor Wheels
- DC TiMotion Motor System
- Side rail and battery backup as option
- Steel Frame with Special Anti-Bacterial Coating



Product Specifications

Max Bed Height	70cm
Min Bed Height	30cm
Length	219cm
Width	121.7cm
Trendelenburg & Rev	-10° - +10°
Capacity	220kg
Product Weight	96kg

Product Specifications

Surface Length	198cm
Surface Width	105cm
Backrest Adjustment	0 - 73°
Knee- Break Adjustment	0 - 35°

ARTG Listing No.

- 196516