

Terry Cloth Sock Aid



Pulling tabs with loops for one-handed application

Terry cloth grips the sock; nylon lining slides up the foot

A terry cloth covered sock aid to help easily apply socks without the need to bend.

The flexible design of the plastic core with slits helps the product to contour around the heel.

Inner nylon lining helps reduce friction while the outer terry cloth lining helps to hold the sock firmly in place until on the foot. Two pulling tabs with loop handles for easy, one-handed application.

- Designed to fit any size sock.
- Machine washable.

Code	PAT-A7273
Length	230 mm
Weight	73 g
Warranty	1 year

